

CareSens™

CareSens System Convenience

Results as fast as five seconds in a unique easy-to-use system

Two test strip ports (top and bottom)

- **Test Range:** 1.1 – 33.3 mmol/L
- **Sample Size:** 0.5 ul
- **Test Time:** 5 seconds
- **Test Sample:** Fresh capillary whole blood
- **Calibration:** Plasma-equivalent
- **Assay Method:** Electrochemical
- **Memory:** 250 blood glucose values

Operating Range

- **Temperature:** 10-40°C
- **Relative Humidity:** 10-90%
- **Hematocrit:** 20-60%

Lifetime
Warranty



CareLance

Vacuum assisted auto-lancing device

Clean, safe, painless retrieval at
your choice of sample site

CareLance
sold separately

Alternative site monitoring with vacuum assisted sampling



Arm



Forearm



Palm



Thigh

fast • painless • accurate

FREE meter to pc download software

