

Result Interpretation (Cat)

Diagnosis of Feline hyperthyroidism and treatment monitoring

< 0.8 µg/dL (< 10.30 nmol/L)	0.8~4.7 µg/dL (10.30 – 60.49 nmol/L)	2.3~4.7 µg/dL (29.60 – 60.49 nmol/L)	> 4.7 µg/dL (> 60.49 nmol/L)
Low	Normal	Gray zone	Consistent with hyperthyroidism

* 1 µg/dL is equal to 12.87 nmol/L



Diagnosis
& Monitoring

T4 Test

< 0.8 µg/dL
(< 10.30 nmol/L)

Euthyroid sick or iatrogenic

0.8 - 2.3 µg/dL
(10.30 – 29.60 nmol/L)

Hyperthyroidism unlikely

2.3 - 4.7 µg/dL
(29.60 – 60.49 nmol/L)

Additional tests including free T4 suggested

> 4.7 µg/dL
(> 60.49 nmol/L)

Hyperthyroidism likely