


Vcheck M

Feline Respiratory 7 Panel


Specification



Detection Target : ① Feline herpesvirus type 1
 ② Feline calicivirus
 ③ Influenza virus
 ④ *Bordetella bronchiseptica*
 ⑤ *Mycoplasma felis*
 ⑥ *Chlamydia felis*
 ⑦ *Streptococcus zooepidemicus*

Species : Feline
Sample Type : Conjunctival swab and Deep pharyngeal swab
Testing Time : 70 minutes
Shelf Life : 12 months

Test Procedure



Performance

Analytical sensitivity (Limit of detection, LoD)

	FHV-1	FCV	Influenza virus	<i>B. bronchiseptica</i>	<i>M. felis</i>	<i>C. felis</i>	SEZ
LOD (copies/ml)	3.3x10 ²	3.3x10 ²	6.7x10 ²	3.3x10 ²	6.7x10 ²	3.3x10 ²	3.3x10 ²

Ordering Information

Product No.	Product Name	Operating Condition	Packing Unit
VCM11AC	Vcheck M Feline Respiratory 7 Panel	2 ~ 28 °C	5 Tests/Kit